

Mass and Hanger Set


ME-8979

Included Parts	Quantity
100 g mass	3
50 g mass	3
20 g mass	6
10 g mass	3
5 g mass	3
2 g mass	3
1 g mass	3
0.5 g mass	3
Mass Hanger	4
Storage Box	1

Replacement Kits	Model Number
Spares Kit including: 3 × 2 g masses 3 × 1g masses 3 × 0.5 g masses 4 × mass hangers	ME-8980
5 g masses, set of 5	ME-8981
10 g masses, set of 5	ME-8982
20 g masses, set of 5	ME-8983
50 g masses, set of 5	ME-8984
100 g masses, set of 5	ME-8985

String Attachment

Pull string into notch and wrap a few times.



Technical Support

Address: PASCO scientific
10101 Foothills Blvd.
Roseville, CA 95747-7100

Phone: 916-786-3800 (worldwide)
800-772-8700 (U.S.)

Fax: (916) 786-7565

Web: www.pasco.com

Email: support@pasco.com

Limited Warranty For a description of the product warranty, see the PASCO catalog.

The PASCO scientific 012-10141A *Mass and Hanger Set Instruction Sheet* is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

PASCO and PASCO scientific are trademarks or registered trademarks of PASCO scientific, in the United States and/or in other countries.